



Fully programmable constant-current electrical stimulator for functional electrical stimulation

- · Delivers biphasic constant current pulses
- Stand-alone device that is configured and controlled from a computer system
- Can send triggers to other devices for synchronization
- Includes electrode impedance check and stimulation current monitoring
- CE-certified (Europe) and FDA-cleared (USA) medical device for use in humans





Authorized Distributor: Ultra Technology Engineering

1115, Block A4, Pusat Dagang Setia Jaya, 9 Jalan PJS 8/9, Bandar Sunway 46150 Selangor, Malaysia Tel:

+603-8066 6726 Email: enquiry@ute-my.com Web: www.ute-my.com









g. Estim FES description

The g.Estim FES is a constant-current, biphasic electrical stimulator intended for powered muscle stimulation, g.Estim FES has an applied part of type BF with connectors for bipolar stimulation electrodes (anode and cathode). The device is controlled by a computer via USB connection. It also has digital outputs for synchronization with other devices. An optional hand-switch allows you to perform stimulation manually. Alternatively, an optional footswitch can be used to explicitely enable/disable stimulation. g.Estim FES includes an impedance measurement function and contains a high-impedance stop mechanism that stops stimulation when electrode contact becomes loose. The device is CE certified (Europe) and cleared by the FDA (USA) for therapeutic use in human patients.

The intended use

The g.Estim FES is a powered muscle stimulator and intended for

- · temporary relaxation of muscle spasms,
- · prevention or retardation of disuse atrophy,
- · increasing local blood flow in the treatment area,
- · muscle re-education,
- · prevention of post-surgical phlebo-thrombosis through immediate stimulation of calf muscles.
- · maintaining or increasing range of motion.

The device must be used by trained and qualified personnel.

General specifications

Stimulus current output Phase shape Work time

Rest time Phase duration

Pulse rate

Rising/falling ramps Power supply

± 1 - 60 mA (±10%) rectangular

2 - 20 s in 0.5 s increments (± 10 % accuracy) 2 -

50 s in 0.5 s increments (± 10 % accuracy) 50 µs - 400 µs in 10 µs increments

 $(\pm 10\% \text{ or } \pm 10 \mu \text{s whichever greater})$

1 - 100 pulses/second in 1 pulse increments

 $(\pm 10\%)$

(Pulse period from 1 s down to 10 ms)

0 - 10 s in 0.5 s increments (± 5 % accuracy) 2 x 9

V battery, USB - connection

The device is designed and manufactured according to the following norms: IEC 60601-1, IEC 60601-1-2, IEC 60601-2-10, IEC 62304, IEC 62366, ISO 14971

Rx only



rear view



